EBA 200

Compact clinical centrifuge



features

Run-Silent™ Operation

Excessive noise in the laboratory can lead to distractions, stress, anxiety, high blood pressure and chronic fatigue. Our centrifuges are as quiet as a normal conversation and will not increase the noise levels within your lab.

© Easy-Lift™

Our Easy-Lift Lid, featuring stay-in-place positioning and lid drop protection, is ergonomically designed to limit repetitive stress injuries, spillage and accidental sample disruption when loading and unloading your centrifuge.

Our Safety-Lock™
Our Safety-Lock Lids are self-closing, auto-locking and come with either single and dual steel lid locks. This ensures a safe, tight closure of the centrifuge chamber during operation. The centrifuge will not operate if lid is not secure. For added safety, our centrifuges feature all-metal containment and adhere to all national and international safety standards and certifications.



overview

The EBA 200 high-performance clinical centrifuge, offers simplicity and safety within a compact, elegant frame.

For convenient operation, the EBA 200 features a built-in, high-performance, 8-place, fixed-angle rotor and a bright, updated LCD display. The EBA 200 is ideal for small-batch clinical chemistry applications.

techspecs

Power source 208 - 240V (50-60Hz) / 100 - 127 V (50-60Hz)

Max. capacity 8 x 15 mL

Max. RPM (speed)

6000 min⁻¹

Max. RCF

3461

Dimensions (H x W x D)

228 x 262 x 352 mm

Weight

approx. 9kg

Ambient consitions

temp: 2°C to 40°C

humidity: max 80% at 31°C - 50% at 40°C

rotorsheet

See our separate EB200 rotor & adaptor sheet for a overview of important tubes.

