## Alaris GP Plus



# Volumetric pump





### features

- Large Display Quickly read critical appropriate data from the large clear display with fewer key presses.
- Colour Coded Set Loading Safe and intuitive set-up with colour coded IV infusion set loading.
- **Door Removal -** With the use of a tool, the door can be safely removed in seconds for cleaning and maintenance.
- Air-in-line Sensor Air-in-line sensor provides both single bubble and accumulated bubble protection.
- Alarm & Warning Beacon An immediate identification of the infusion pump requiring attention.
- User Friendly Interface Featuring the same intuitive user interface as the Alaris™ Syringe Pumps, the GP Plus volumetric pump requires minimal additional training.

#### overview

Acute and sub-acute applications require a simple and easy to use volumetric pump. Our Alaris® GP volumetric pump provides a versatile and flexible platform for infusion therapy which is suitable across a broad range of applications.

The Alaris® GP Volumetric pump offers a range of features suited to drug therapy, blood transfusions and parenteral feeding: including a large, clear specialty sets.

### techspecs

#### Infusion Rate Range

0.1 - 1200ml/h (with titration if enabled) **VTBI** 

Range: Primary - 1ml to 9999ml End Actions: Stop, KVO and continue

Infusion Rate Range

0.1 - 1200ml/h (with titration if enabled)

Volumetric Accuracy +/- 5% (nominal)

Hands on and Hands Free Bolus

Bolus Rate 10 - 1200 ml/h Bolus Volume 0.1 - 100 ml

Display

Mono (240x240 pixels)

**KVO Rates** 

0.1-20ml/h or set infusion rate if lower



